


Mammography



Extras

Any Brand

Not brand specific - works through DICOM

Workstation

Computer for image analysis, also called postprocessing. Has two screens of 5mp each to reach 50 micron resolution per pixel. The system is functional without a workstation, but it is highly recommended for cancer screening.



Printer

A high definition printer that is adequate for mammograms. 50 microns is the guideline. Only necessary if you need pictures in print.



Brand Specific

Has to be same brand as system

Tomosynthesis

Also called 3D, is an option that can be added. However, it requires that the system, from factory, has the hardware for 3D.



Stereotactic

Biopsy kits are an extra option that eases the process of obtaining samples in case of tissue changes identified on the mammogram.

